

LOQs for Damask Rose by LC-MS/MS

No.	Compound Name	LOQ (ppb)
1	Abamectin	10
2	Acephate	10
3	Acetamiprid	10
4	Acibenzolar-s-methyl	10
5	Aldicarb Sulfone	10
6	Ametryn	10
7	Aminocarb	10
8	Azinphos-methyl	10
9	Azoxystrobin	10
10	Benalaxyl	10
11	Bendiocarb	10
12	Benzoximate	10
13	Bifenazate	10
14	Buprimate	10
15	Butoxycarboxim	10
16	Carbetamide	10
17	Carbofuran	10
18	Carbofuran 3-Hydroxy	10
19	Carboxin	10
20	Carfentrazone-ethyl	10
21	Chlorantraniliprole	10
22	Chlorbromuron	10
23	Chlorfluazuron	10
24	Chloroxuron	10
25	Chlortoluron	10
26	Clodinafop-propargyl	10
27	Clofentezine	10
28	Clothianidine	10
29	Coumachlor	10
30	Coumafuryl	10
31	Coumaphos	10
32	Coumatetralyl	10
33	Crimidine	10
34	Cyanofenphos	10
35	Cyazofamid	10
36	Cymoxanil	10
37	Cyproconazole I,II (Total)	10
38	Cyprodinil	10
39	Desmedipham	10
40	Diclobutrazole	10
41	Diclofluanid	10
42	Dicrotophos	10
43	Diethofencarb	10
44	Diflubenzuron	10
45	Dimethachlor	10
46	Dimethoate	10
47	Dimethomorph I, II (Total)	10
48	Dimoxystrobin	10
49	Dinotefuran	10
50	Dioxacarb	10
51	Disulfoton sulfone	10
52	Disulfoton sulfoxide	10
53	Doramectin	10
54	Epoxiconazole	10

No.	Compound Name	LOQ (ppb)
55	Etaconazole	10
56	Ethiofencarb	10
57	Ethofumesate	10
58	Ethoprophos	10
59	Fenamidone	10
60	Fenamiphos	10
61	Fenbucarb	10
62	Fenbuconazole	10
63	Fenoxycarb	10
64	Fenoxyp-ethyl	10
65	Fenpropimorph	10
66	Fenpyroximate	10
67	Fenuron	10
68	Fonicamid	10
69	Fludioxinil	10
70	Flufenacet	10
71	Flufenoxuron	10
72	Fluometuron	10
73	Fluquinconazole	10
74	Flusilazole	10
75	Flutolanil	10
76	Forchlorfenuron	10
77	Formothion	10
78	Fuberidazole	10
79	Furathiocarb	10
80	Halofenozide	10
81	Haloxifop-methyl	10
82	Hexaflumuron	10
83	Hexazinon	10
84	Hexythiazox	10
85	Hydramethylnon	10
86	Imazalil	10
87	Imidacloprid	10
88	Indoxacarb	10
89	Iaconazole	10
90	Iprovalicarb	10
91	Isoprocarb	10
92	Linuron	10
93	Lufenuron	10
94	Mandipropamid	10
95	Mecarbam	10
96	Mefenacet	10
97	Mepanipyram	10
98	Mepronil	10
99	Metaflumizone	10
100	Metconazole	10
101	Methabenzthiazuron	10
102	Methamidophos	10
103	Methiocarb	10
104	Methomyl	10
105	Methoprotryne	10
106	Metobromuron	10
107	Mevinphos I, II (Total)	10
108	Mexacarbate	10

No.	Compound Name	LOQ (ppb)
109	Monocrotophos	10
110	Monolinuron	10
111	Moxidectin	10
112	Neburon	10
113	Novaluron	10
114	Nuarimol	10
115	Omethoate	10
116	Oxadixyl	10
117	Oxamyl	10
118	Oxydemeton methyl	10
119	Paclbutrazole	10
120	Pencycuron (monceren)	10
121	Pendimethalin	10
122	Phenmedipham	10
123	Phosmet	10
124	Picoxystrobin	10
125	Prochloraz	10
126	Promecarb	10
127	Prometon	10
128	Prometryn	10
129	Propham	10
130	Propoxur	10
131	Pyracarbolid	10
132	Pyraclostrobin	10
133	Pyrimethanil	10
134	Quinoxifen	10
135	Siduron	10
136	Simetryn	10
137	Spinetoram	10
138	Spiromesifen	10
139	Spirotetramate	10
140	Spiroxamine I, II (Total)	10
141	Sulfaquinoxaline	10
142	Sulfentrazone	10
143	TCMTB	10
144	Tebufenpyrad	10
145	Tebufenzoid	10
146	Tebuthiuron	10
147	Temephos	10
148	Terbumeton	10
149	Tetraconazole	10
150	Thiacloprid	10
151	Thiamethoxam	10
152	Thiobencarb	10
153	Thiodicarb	10
154	Triazophos	10
155	Tridemorph	10
156	Triflumizole	10
157	Triflumuron	10
158	Triticonazole	10
159	Uniconazole	10
160	Vamidathion	10
161	Warfarin	10
162	Zoxamide	10